Dentsply Sirona

Enlargers

EN	Instructions For Use	2
FR	Instructions d'utilisation	7
IT	Istruzioni per l'uso	12
DE	Gebrauchsanweisung	17
ES	Instrucciones de uso	22
PT	Instruções de utilização	27
NL	Gebruiksinstructies	32
SV	Användarinstruktioner	37
DA	Brugsanvisning	42
NO	Brukerinstruksjoner	47
LT	Naudojimo instrukcija	52
LV	Lietošanas norādījumi	57
ET	Kasutusjuhised	62
CS	Návod k použití	67
SK	Návod na použitie	72
HU	Használati utasítás	77
RO	Instrucţiuni de utilizare	82
SR	Uputstva za upotrebu	87
BS	Uputstva za upotrebu	92
CNR	Uputstva za upotrebu	97
SQ	Udhëzimet për përdorim	102
EL	Οδηγίες χρήσης	107
TR	Kullanım Kılavuzu	112
PT-BR	Instruções de uso	116
KO	사용 지침	121

Dentsply Sirona

EN

Enlargers

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE
NON-STERILE – MANDATORY PROCESSING BEFORE FIRST USE
SINGLE USE – REPROCESSING PROHIBITED
REUSABLE – LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE NON-STERILE ROTARY INSTRUMENTS
REF. A0007XXXXXX00 (Root Canal Drill), A0008XXXXXX00 (Gates Glidden Drill, X-Gates), A0009XXXXXX00 (Largo® Peeso Reamer)
SINGLE USE STERILE ROTARY INSTRUMENTS
REF. A0008XXXXXX03 (Gates Sterile Glidden Drill, X-Gates Sterile), A0009XXXXXX03 (Largo® Sterile Peeso Reamer)

1 PRODUCT DESCRIPTION

- Product description:
 Enlargers are single piece instruments made of stainless steel with three different parts:
 - The active part,
 - The neck,
 - The 13 mm shank with notch (for instrument size identification).
- Clinical benefits:
 The clinical benefit for the patient being treated with Enlargers is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Enlargers are engine driven instruments intended for root canal preparation (enlargement of the coronal part / preflare / preparation of the entrance).
- Intended patient population:
 Patients who need root canal treatment.

- Intended user:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

Enlargers can be delivered in different versions: sterile blisters or non-sterile plastic boxes (cassettes).

The available types and sizes of Enlargers instruments are detailed in the following table:

SKU	Product	Available lengths	Available size N°X Ø (number of notches)	Packaging and content
A0007	Root Canal Drill	33 mm	N°1 Ø 1.20 mm (1 notch) N°2 Ø 1.40 mm (2 notches) N°3 Ø 1.65 mm (3 notches) N°4 Ø 1.90 mm (4 notches)	Box of 12 cassettes: each including 6 N°X instruments
A0008*****X00	Gates Glidden Drill	28 mm 32 mm 34 mm	N°1 Ø 0.50 mm (1 notch) N°2 Ø 0.70 mm (2 notches) N°3 Ø 0.90 mm (3 notches) N°4 Ø 1.10 mm (4 notches) N°5 Ø 1.30 mm (5 notches) N°6 Ø 1.50 mm (6 notches)	Box of 12 cassettes: each including 6 N°X instruments
A0008*****X03	Gates Sterile Glidden Drill	28 mm 32 mm		Box of 8 blisters: each including 6 N°X instruments
A000825090100	Gates Glidden Drill	34 mm		Box of 12 cassettes, each including 6 instruments: 1 N°1 + 1 N°2 + 1 N°3 + 1 N°4 + 1 N°5 + 1 N°6
A000824090103	Gates Sterile Glidden Drill	28 mm 32 mm		Box of 8 blisters, each including 6 instruments: 1 N°1 + 1 N°2 + 1 N°3 + 1 N°4 + 1 N°5 + 1 N°6
A0008*****700	X-Gates	32 mm	Ø 1.10 mm	Box of 12 cassettes: each including 6 X-Gate instruments
A0008*****703	X-Gates Sterile			Box of 8 blisters: each including 6 X-Gate Sterile instruments
A0009*****X00	Largo® Peeso Reamer	28 mm 32 mm	N°1 Ø 0.70 mm (1 notch) N°2 Ø 0.90 mm (2 notches) N°3 Ø 1.10 mm (3 notches) N°4 Ø 1.30 mm (4 notches) N°5 Ø 1.50 mm (5 notches) N°6 Ø 1.70 mm (6 notches)	Box of 12 cassettes: each including 6 N°X instruments
A0009*****X03	Largo® Sterile Peeso Reamer			Box of 8 blisters: each including 6 N°X instruments
A000924090100	Largo® Peeso Reamer			Box of 12 cassettes, each including 6 instruments: 1 N°1 + 1 N°2 + 1 N°3 + 1 N°4 + 1 N°5 + 1 N°6
A000924090103	Largo® Sterile Peeso Reamer			Box of 8 blisters, each including 6 instruments: 1 N°1 + 1 N°2 + 1 N°3 + 1 N°4 + 1 N°5 + 1 N°6

5 COMPOSITION

Enlargers are single piece instruments made of Stainless steel.

6 COMPATIBLE DEVICES

Enlargers shall be connected to contra-angles and used in combination with endodontic motors at constant rotation speeds as follow:

- X-Gates: within the 500 - 800 rpm range,
- Gates Glidden Drill: within the 800 - 1200 rpm range,
- Root Canal Drill: within the 800 - 1200 rpm range,
- Largo® Peeso Reamer: within the 800 - 1200 rpm range.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the devices, the patient and/or the user:
 - Serious health consequences (infection, toxic or allergic reactions)
 - Swallowing or inhalation of the instrument
 - Tooth extraction
 - Tissue necrosis
 - Discomfort
 - Treatment delay
- The use of inappropriate rotation speed, pressure and/or contact duration would speed up the wear of the instruments, increase the risk of breakage and heat creation.
- Sterile instruments:
 - Sterility cannot be guaranteed if the packaging is open, damaged or wet.
 - Do not use sterile devices after expiration date.
 - Single use devices cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Non-sterile instruments:
 - Processing step is mandatory before use for non-sterile devices mentioned in this IFU (see [8]).

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask, waterproof gown).
- Inspect the packaging before use, do not use and discard the devices if the packaging is damaged.
- Inspect the device before each use for signs of defects such as dull spots, deformations (bent), breakage, corrosion or damaged cutting edges. With these indications, the device is not able to fulfil the intended use with the required safety level and should be discarded.
- Before using any device, make sure it is well connected to the handpiece.
- Do not apply excessive pressure on the device.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.

9 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

10 EXPIRATION DATE

Do not use sterile devices after the expiration date indicated on the label.

11 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

12 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

13 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE
2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH